Investigating Leaky Gut Testing the Gut Barrier

Emily Birch
KBMO Diagnostics



What are we covering today?

- Overview of 'leaky gut'
- How to test the gut barrier
- The food sensitivity gut barrier connection
- Gut barrier markers
- Real-life client reports



Meet the Team



Charlotte Hunter Head of KBMO UK



Clare Kennedy Operations Manager



Linette Petrillo Customer Services



Natasha Khan Practitioner Relationships



Kelly Hutson Business Relationships




Emily Birch Support & Communications






"All disease begins in the gut" - Hippocrates



Source: Afzaal et al. 2022



The Intestinal Barrier (IB)




Source: Sharma & Riva, 2022

- The first line of defence against pathogens and food antigens
- In a healthy person, it operates on selective permeability. Water, ions and low-molecular substances can freely pass through.
 Macromolecules, toxins, food allergens and pathogens cannot.
- When the latter substances permeate the IB, the immune system is over-activated, resulting in a systemic inflammatory response.
- Maintaining the integrity of the intestinal barrier is key for your clients' wellbeing.



Intestinal Permeability – "Leaky Gut"



Source: Kocot, Jarocka-Cyrta & Drabinska, 2022

- The lining of the small intestine becomes damaged and more porous.
- Allows toxins, bacteria and other particles to "leak" into the bloodstream.
- Contributes to excessive immune reactions.
- Triggers and mediates systemic inflammation.
- Is associated with a wide range of symptoms, including food sensitivities, bloating, fatigue, headaches, joint pain, skin issues and nutritional deficiencies.

Remember: A level of intestinal <u>cell</u> permeability is normal for nutrient absorption! But "leaky gut" involves compromised epithelial tight junctions.



Testing the Gut Barrier



The Food Inflammation Test (FIT)






Dr Joel Evans

"The immune complexes and inflammation associated with foods are an underlying problem for all of the conditions I see daily."




Food sensitivity or leaky gut: which came first?





Why test food sensitivities?





I know what you're probably thinking...





Why is isolated IgG testing problematic?




High rate of false positives





Technology for Enhanced Sensitivity





The Immune System

The immune system provides 3 levels of defence against disease causing organisms.

The FIT and Gut Barrier Panel tests all of these:

1. Barriers

Prevent Entry

- · Skin & mucus membranes
- · Stomach acid & digestive enzymes
- Beneficial bacteria that live in the colon (the gut microbiota)

2. Innate Immunnity

General Defence

 WBCs called neutrophils & macrophages engulf & destroy foreign invaders & damaged cells

3. Adaptive Immunity

Specific Defence

- WBCs (T lymphocytes or T cells) target & destroy infected or cancerous cells
- WBCs called B lymphocytes (B cells) & plasma cells produce antibodies that target & destroy infected /cancerous cells



Foods We Test

			176 Foods & Gut Ba	rrier Panel			
132 Foods & Gut Barrier Panel							
		:	22 Foods				
DAIRY/EGGS Casein Cow's Milk Egg White Egg Yolk MEATS Beef Chicken	FISH Salmon SEAFOOD Shrimp	GRAINS Wheat, Gluten Wheat, Whole BEANS Coffee Soy Bean	VEGETABLES Com Potato, White Tomato	FRUITS Banana Pineapple	NUTS/SEEDS Almond Peanut	SPICES/ MISCELLAN Candida Turmeric Yeast, Brewer's	NEOUS
Goat's Milk Whey Bacon Duck Lamb Pork Turkey	Codfish Flounder Halibut Sea Bass Snapper Swordfish Trout Tuna Clam Crab Lobster Scallops	Barley Millet Oat Quinoa Rice Rye Cocoa Kidney Bean Lentils Navy Bean Pinto Bean	Artichoke Asparagus Pumpkin Beets Spinach Broccoli Zucchini Butternut Squash Cabbage Carob Carrot Cauliflower Celery Collard Greens Cucumber Lettuce Pea, Chick Pea, Green Pepper, Green	Apple Avocado Blueberry Cantaloupe Cherry Cranberry Grape, White Grapefruit Honeydew Melon Lemon Lime Olive, Green Onion, White Orange Peach Pear	Cashew Coconut Cola Nut Dill Seed English Walnut Flax Seed Hazelnut Pecan Sesame Seed Sunflower Seed	Agave Aspartame Basil Benzoic Acid BHA Canola Oil Cinnamon Garlic Ginger Hops MSG Mushroom Mustard Oregano Paprika Pepper, Black	Pepper, Chili Peppermint Polysorbate 80 Red #3 Red #40 Rosemary Saccharin Spirulina Sugarcane Tea, Black Vanilla Wine, Red Yellow #6
Sheep's Milk Venison	Anchovy Mackerel Sardine Oyster Squid	Amaranth Buckwheat Gliadin Sorghum Spelt Black Bean Green Bean	Arugula Brussel Sprouts Cilantro Kale Parsley Summer Squash	Acai Berry Apricot Eggplant Fig Kiwi Goji Berry Mango Monk Fruit Papaya	Brazil Nut Chia Seed Hemp Seed Macadamia Nut Pine Nut Pistachio	Cloves Coconut Oil Cumin Hemp Protein (O Honey Maple Syrup Stevia Tapioca Vinegar	CDB)

Food selection contains 'real-world' foods, encompassing raw and cooked.

All FIT tests include the Gut Barrier Panel.



Summary of FIT Tests

FIT22 22 of the most common food sensitivities including gluten, cow's milk and egg plus the Gut Barrier Panel

FIT132

132 foods and food additives plus the Gut Barrier Panel


FIT176


176 foods including health foods such as honey, stevia and coconut oil plus the Gut Barrier Panel

Assess the integrity of the gut lining by measuring Candida, Zonulin, Occludin and Lipopolysaccharides (LPS)



The FIT Report





- Easy to read
- Colour coded
- Client-friendly
- Easy to interpret
- Efficient practice
- Easy to repeat tests




Easy Interpretation

Wheat, Gluten
Wheat, Whole

Sheep's Milk
Rye
Pomegranate
Carrot
Honey

Mango
Asparagus
Vanilla
Squid



Each food is colour-coded and graded based on reactivity intensity:




Red	+4	Highest immune reactivity	
Orange	+3	Strong reactivity	
Yellow	+2	Moderate reactivity	
Light Green	+1	Low grade reactivity	
Dark Green	+0	No reactivity (safe to consume)	



The Gut Barrier Panel



Gut Barrier Panel						
		IgG1-4+C3d	IgA1-2			
		Cut off		Cut off		
Candida	Negative		Positive			
Zonulin	Negative		Positive			
Occludin	Negative		Negative			
LPS	Negative		Positive			



Our unique Gut Barrier Panel acknowledges that leaky gut occurs across a spectrum and includes the gatekeeper markers: Candida, Zonulin, Occludin and LPS.

For each marker we measure IgG 1-4 / C3d, in addition to IgA 1-2.



Candida

Occludin

Naturally occurring yeast, residing in the GI tract as part of the normal microbiome. Overgrowth is problematic.

We measure and use any candida overgrowth in the stomach/dysbiosis as a precursor to leaky gut occurrence

A marker of intestinal tight junction stabilization and optimal barrier function.

Elevated occludin indicates that the tight junctions between intestinal epithelial cells are breaking down

Zonulin

A marker of intestinal permeability. We use unique next generation Zonulin IgG **antibody** screening – a more stable and specific marker, exclusive to KBMO.

Developed by Dr Alessio Fasano & Dr Brent Dorval

Lipopolysaccharide (LPS)

Potent endotoxin, present in outer surface membrane of gram-negative bacteria - many of which are pathogenic.

Major inducer of inflammatory response triggers inflammatory cytokine release, can create direct epithelial damage in the gut, crosses the blood brain barrier.



Real Clients, Real Reports But before we dive in...

How strong is your rapport with your client?

Will you need to negotiate with your client?

Are there any clientspecific contraindications? How will you decide what to do with the results?

How will you coach your client?

What's the clinical context?

DO NOT look for things that give you an excuse to treat!

Your client is unique any 'protocols' require personalisation



Client 1 – Anna

Female, 47 years, post-menopausal (surgical menopause – adenomyosis & endometriosis)


History of: Anxiety, controlling behaviour around foods

Symptoms: Postprandial bloating, constipation, reflux, brain fog, generalized fatigue, joint pains, frequent headaches

Prior testing: SIBO test (negative), H.Pylori test (negative), GI360 Stool Analysis (unremarkable), various blood tests (all within normal range), Coeliac testing (normal)


Diet: Reasonably healthy, no alcohol, good hydration, balanced macros, little sugar

Lifestyle: Daily exercise, low reported stress





Client 1 - Anna







Gut Barrier Panel							
		lgG1-4+C3d	IgA1-2				
		Cut off		Cut off			
Candida	Negative		Negative	<u> </u>			
Zonulin	Negative		Positive				
Occludin	Negative	<u> </u>	Negative	<u> </u>			
LPS	Negative		Negative				

Longer term, systemic


Short term, immune activation

C3d:

- Downstream inflammatory metabolite
- Breakdown product of complement component C3
- Can indicate complement pathway activation in response to immune complexes.



Client 1 - Anna



What we did:

Food: Exclusions as per her FIT Test results. Healthy replacement of gluten. Increased polyphenols and fibre (with as much variety as possible).

Lifestyle: Addressed eating 'housekeeping' (sufficient chewing, limited water at meals, mindful eating), increased sleep quantity

Supplements: Nutridyn: GI Integrity Peach Tea, Nutridyn: Fiber

Don't assume candida is the problem, just because its there.

Be careful with your hand grenades!




Client 2 - Christine

Female, 69 years, prior breast cancer history

Symptoms: Severe allergies, histamine intolerance, fatigue, anxiety, Celiac disease





Approach: Continued gluten & dairy-free diet, gut healing (digestive enzymes, Mg, Zn, B-vitamins), gentle anti-Candida protocol, adaptogens and adrenal cortex glandulars.

Progress at Follow Up: All GBP markers negative. Improved energy, digestion and allergies. HbA1c improved from 5.8% to 5.1%.



Client 3 - Mark



Male, 30 years

Symptoms: Ulcerative Collitis flare, loose stools with urgency, fatigue, poor exercise performance, increased CRP

	Spring 2024		11 months		End of	Winter 2026	=
	Spring 2024 —			End of Winter 2025			
	Gut Barrier Panel				Gut Barrier Pane	el	
	IgG1-4+C3d	IgA1-2			lgG1-4+C3d		IgA1-2
Candida	Positive	Positive	- "	S. Charles	Cut off	Name of the last	Cu
Zanulia	Negative	Positive	Candida	Negative		Negative	
Zonulin	Negative	Positive	Zonulin	Negative		Negative	
Occludin	Positive	Negative	Occludin	Negative		Negative	
LPS	Positive	Positive	LPS	Negative		Negative	

Approach: Removed gluten, dairy and FIT-identified foods. Added whole foods (++ polyphenols). Gut barrier healing protocol, betaine HCl and digestive enzymes, anti-inflammatory herbs.

Progress at Follow Up: Symptoms resolved. Energy strong. Digestion normal. HbA1c improved. UC flare fully resolved without steroids (client's choice).



Gut Barrier Support



Optimise dietary fibre - feed the microbes so they don't eat your mucins!



Polyphenol-rich diet - berries, pomegranate, blackcurrant, green tea, cacao



Gut Supporting nutrients - Glutamine, Glycine, Zinc, EFA's, B vitamins



Mucopolysaccharides / glucosaminoglycans - Bone broth, fish/poultry skin, flaxseeds, aloe vera, oats, okra, medicinal mushrooms, apples, citrus peel, slippery elm, marshmallow



Address food sensitivities / food-mediated inflammation




Mindful eating




Stress reduction, sleep optimisation, breathing techniques



Supplement Guidance









When to re-test

We recommend re-testing after 3 – 6 months




What if I don't see changes, or they get worse?

- See previous FIT test considerations
- Consider ACTUAL client compliance
- Consider where you are in your support protocol
- Consider increased stress and other lifestyle or health factors
- Consider running a stool / SIBO test
- Consider GP referral or involvement



Clinical Support



- East to read patient reports
- Educational videos and webinars
- Provider handbook
- Supplementary handouts and leaflets
- Supplement protocols
- Website resources
- Support calls













Patient Support



- East to read patient reports
- Client Mobile App
- Meal suggestions





A powerful tool for your practice

- Useful pre- initial consultation screening tool for every client
- Data you can use straight away (immediate clarity)
- Improves client engagement & compliance
- No guessing at problematic foods / gut barrier compromise
- It's not a "so what?" test, it's insightful and useful
- Financial benefit for you, with a practitioner mark-up of your choice
- Clear, concise reports mean less need for long support calls and assistance = freeing up your time

The gut is the foundation of systemic health

It influences immune regulation, detoxification, hormone balance, and brain function.

The KBMO Gut Barrier Panel can help you identify hidden dysfunction, including silent immune activation





"You cannot repair what you haven't measured. Gut barrier testing gives us that window – and the data to guide true healing."

- HUM2N Clinical Team



Stay in Touch



ukoffice@kbmodiagnostics.com



www.kbmodiagnostics.co.uk



@kbmodiagnosticsuk





Thank You



Image References

Afzaal, M. Saeed, F. Shah, Y.A *et al.* (2022). 'Human gut microbiota in health and disease: Unveiling the relationship', *Frontiers in Microbiology,*

Kocot, A.M. Jarocka-Cyrta, E. & Drabinska, N. (2022). 'Overview of the Importance of Biotics in Gut Barrier Integrity', *International Journal of Molecular Sciences*, 23 (5).

Sharma, L. & Riva, A. (2022). 'Intestinal Barrier Function in Health and Disease – Any Role of SARS-CoV-2?', *Microorganisms*, 8 (11).

